



# The Static SIMS Library Version 4

Secondary Ion Mass Spectraデータベース



SIMSの(Secondary Ion Mass Spectraデータベースライブラリー。約1000物質以上についてPositive, Negative ion spectraを約1900スペクトル収録しました。製作会社より提供されたスペクトル・データベースは標準フォーマットに変換されオリジナル・スペクトルの詳細すべてがオリジナルのデータに使われていたマス・レゾリューションでアクセスすることができます。

Material	CAS Registry No.	Meckl Index No.	Mol. Wt.	Particle	Sample conditions
[16] Dendrimer			1676.36	Cs+	Spin cast from chloroform onto glass
[19] Dendrimer			2238.43	Cs+	Spin cast from chloroform onto glass
[7] Dendrimer			798.17	Cs+	Spin cast from chloroform onto glass
[L]ucine 5[en]kephalin (cationised)	58822-25-6	11.3528	555.27	Ar+	2 [E]l of solution was deposited on t
[M]ethionine 5[en]kephalin (cationised)	58569-55-4	11.3528	573.23	Ar+	4 [E]l of solution was deposited on t
1,4-Dihydroxyethylol	6892-63-8		154.01	Ar+	Powder pressed into indium
12-Hydroxyoctanoic acid	106-14-9		300.27	69 Ga+	Powder
15-Mercapto-1-carboxymethylpentadecanoic acid				Cs+	SAM on gold
15-Mercapto-1-pentadecanoic acid				Cs+	SAM on gold
16-Mercapto-1-hexadecanoic acid				Cs+	SAM on gold
16-Mercapto-1-hexadecanol				Cs+	SAM on gold
272.14			272.14	69 Ga+	~0.2 nm thick disc. Pressed from
202.18			202.18	Ar+	Chemisorption of mercaptan on sil
448.05			448.05	69 Ga+	Neat film frozen in-situ
350.11			350.11	69 Ga+	Thin film, spin cast from methanol c
382.16			382.16	Cs+	Sample mounted under a metal grid



ライブラリー表示

**データ表示**

Search Criteria

Mass	Polarity
* 153	Pos

Add Remove Clear All

Options

Search 4 ions per 50 amu

Ag cationised materials only

Search Close

**検索画面**

Results

Hits	Material	Mol. Wt.
1	Poly(BA - co - C4), sequential	9.5k
1	Poly[(1,4-bisacryloylpiperazine)(2-methylpiperazine)]	
1	Poly(bisphenol A - co - epichlorohydrin)	
1	Poly(ethylene - co - vinyl acetate) (25% VA)	
1	Poly(ethylene - co - vinyl acetate) (40% VA)	
1	Poly(melamine formaldehyde - co - vinyl alcohol) (20% MF)	
1	Poly(methyl methacrylate - co - PEG 1000 methacrylate), (12.5% PEGMA)	
1	Poly(methyl methacrylate - co - PEG 1000 methacrylate), (50% PEGMA)	
1	Poly[(N,N'-methylene bisacrylamide)(N,N'-dimethyl ethylenediamine)]	
1	Poly(styrene - co - acrylonitrile), (30% AN)	
1	Poly(styrene - co - butadiene - co - styrene), block	
1	Poly(styrene - co - maleic anhydride)	
1	Poly(vinyl alcohol - co - vinyl butyrate)	
1	Poly(vinyl chloride - co - vinyl acetate) (17% VA)	

(Press and hold the Ctrl key to select more than one material)

Load

■ 価格:

新規購入	¥720,000.- (税別)
Ver.2からのアップグレード	¥200,000.- (税別)

【輸入・販売元】

株式会社デジタルデータマネジメント  
〒103-0025 中央区日本橋茅場町1-11-8  
紅萌ビル4F

【販売代理店】

株式会社テクノサイエンス  
〒264-0034 千葉市若葉区原町929-8  
TEL:043-206-0155 FAX:043-206-0188  
URL <https://www.techno-lab.co.jp/>